



## SM-B8

## Base for Stereo Microscope



Model	SM-B8
Base	With clip
	With black & white platen ø95mm
	Size 205x275x42mm
Pillar	Height 240mm
	ø32mm
Reflected illumination	3W double LED lamp, independent switch
	Brightness control
Arm type	SM-A1
Power supply	110/220V 50/60Hz 1pH
Weight	3kg

## SM-B9

## Base for Stereo Microscope



Model	SM-B9
Base	With clip
	With glass platen ø95mm
	Size 300x295x58mm
Draw tube	Height 358mm
	Size ø76mm
Transmitted illumination	3W LED lamp
	Brightness control
Power supply	110/220V 50/60Hz 1pH
Weight	3kg



**DEVCO S.r.l.**

Via Marzabotto, 59 20037 Paderno Dugnano (MI)  
Tel +39 0283591153 - Fax. +39 0295441300  
[www.devcosrl.it](http://www.devcosrl.it) - e-mail: [info@devcosrl.it](mailto:info@devcosrl.it)

[www.echo-lab.it](http://www.echo-lab.it)