

SM-B6

Base for Stereo Microscope



Model	SM-B6
Base	With clip
	With black & white platen ø95mm
	Size 200x255x22mm
Draw tube	Height 300mm
	Size ø76mm
Focus distance	100mm
Weight	1kg

SM-B7

Base for Stereo Microscope



Model	SM-B7
Base	With clip
	With glass platen ø95mm
	Size 205x275x42mm
Drawtube	Height 300mm
	ø76mm
Transmitted illumination	3W LED lamp
	Brightness control
Reflected illumination	3W LED lamp
	Brightness control
Focus distance	100mm
Power supply	110/220V 50/60Hz 1pH
Weight	3kg



DEVCO S.r.l.

Via Marzabotto, 59 20037 Paderno Dugnano (MI)
Tel +39 0283591153 - Fax. +39 0295441300
www.devcosrl.it - e-mail: info@devcosrl.it

www.echo-lab.it