

SM 210 H / SM 216 H /  
SM 215 H / SM 214 H

## SM 210 H

## Stereo Microscope Head




Model	SM 216 H	SM 215 H	SM 214 H
Optical system	Greenough stereoscopic	Greenough stereoscopic	Greenough stereoscopic
Total magnification	20x/40x	10x/30x	10x/20x
Eyepiece	H WF 10x (ø20mm)	H WF 10x (ø20mm)	H WF 10x (ø20mm)
Turret type objective	2x/4x	1x/3x	1x/2x
Zoom mode	2 objectives switchable by eyepiece rotating tube	2 objectives switchable by eyepiece rotating tube	2 objectives switchable by eyepiece rotating tube
Head	Binocular	Binocular	Binocular
	Inclination of 45°	Inclination of 45°	Inclination of 45°
	Interpupillary distance 54-76mm	Interpupillary distance 54-76mm	Interpupillary distance 54-76mm
	Diopter adjustable ±5	Diopter adjustable ±5	Diopter adjustable ±5
Dimensions	310x210x225mm	310x210x225mm	310x210x225mm
Weight	1.5kg	1.5kg	1.5kg

## Optional accessories

Model	Description
SM200eyp015x	Eyepiece WF 15x (ø15mm)
SM200eyp020x	Eyepiece WF 20x (ø10mm)
SM200eyp025x	Eyepiece WF 25x (ø9mm)
SM200210lens0.5x	Objective PL 0.5x WD=165.00mm
SM200210lens1.5x	Objective PL 1.5x WD=45.00mm
SM200210lens2x	Objective PL 2x WD=30.00mm
SM000002	LED ring with 60 bulbs
SM000001	Eyeshield

COMPONENTS DIAGRAM

SM 216 H / SM 215 H / SM 214 H





**DEVCO S.r.l.**

Via Marzabotto, 59 20037 Paderno Dugnano (MI)  
Tel +39 0283591153 - Fax. +39 0295441300  
[www.devcosrl.it](http://www.devcosrl.it) - e-mail: [info@devcosrl.it](mailto:info@devcosrl.it)

[www.echo-lab.it](http://www.echo-lab.it)